


Holds up to heavy expectations

Juno DRF Digital radiography and fluoroscopy designed for bariatrics



PHILIPS
sense and simplicity

When the going gets tough, the tough get going

Philips helps you overcome challenges. One billion overweight adults are a fact we respond to. To conveniently carry out digital radiography and fluoroscopy on patients with obesity, Philips provides you with its Juno DRF 2-in-1 system. A very high weight capacity, the large table-top size and ample distance between table and tube make the system especially suitable for bariatric people. Both patients and medical personnel will experience a more comfortable exam.

A weighty issue

Are you willing to run the risk of sending an overweight patient away because you don't trust the strength of your table? With Juno DRF's 626 lb (284 kg) table capacity this scenario is more than unlikely. A very low table top height and a wide footrest make it easy to position patients with limited mobility. In addition, technologists have space to work around the patient thanks to extended clearance and a high scan range.

Excellent bariatric capabilities support obese patients and ease positioning



A wide range for obesity patients with Juno DRF

- Designated table design with weight capacity of up to 626 lb (284 kg) in all positions combined with large table top
- Low-level table top at 23" (63 cm) facilitates easy transfer and positioning of limited mobility patients
- High clearance with Source to Image Distance (SID) from 43" to 71" (110 cm to 180 cm) provides excellent accessibility
- A scan range of 79.9" (203 cm) without table-top movement enables exams without repositioning
- Large detector size of 17" x 17" (43 cm x 43 cm) ensures optimal patient coverage for different types of exams including lungs and pelvis without repositioning
- The solid footrest supports bariatric patients and makes positioning convenient
- Optional 80 kW generator for shorter exposure times and 800 kHU tube to meet increased tube performance demands

Not available for sale in Canada

Please visit www.philips.com/JunoDRF



© 2010 Koninklijke Philips Electronics N.V.
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Healthcare is part of Royal Philips Electronics

www.philips.com/healthcare
healthcare@philips.com

Printed in The Netherlands
4522 962 66951 * NOV 2010