



Birthday bash or the Bahamas

EverGo has up to 8 hours of battery
life to take you where you want to go

PHILIPS

RESPIRONICS

EverGo

portable oxygen concentrator

Wherever, Whenever

Whether you're off to a birthday party or going on vacation, the EverGo portable oxygen concentrator delivers all-day freedom. Its eight-hour battery life gives you tremendous flexibility to enjoy more of the activities you love. Go wherever you want, whenever you want.

EverGo is ideally suited for people who enjoy travel as it has been approved by the FAA* for use on commercial airlines. Not even the sky is the limit.

* Airline policies for traveling with an oxygen concentrator vary. Patients need to check with their airlines before flying with EverGo.




Live life on your terms



“The quality of life is just fantastic with EverGo. You can do anything. EverGo has just made my life so much easier. I can do whatever I want.”

Dani, actual EverGo user, Michigan



“EverGo definitely helped me. I can do more things. I love it. Couldn't do without it.”

Cliff, actual EverGo user, Pennsylvania



“I would certainly recommend EverGo because it lasts for such a long time. It's wonderful, very easy to use.”

Darlene, actual EverGo user, Florida



“I was pretty confined when I had the tanks. Now, with EverGo, I can go golfing.”

David, actual EverGo user, Ohio

Is EverGo right for me?

EverGo is intended for prescription use by patients requiring oxygen on a supplemental basis. It may be considered for continuous use in home, institutional, and travel/mobile environments. EverGo may not be for all oxygen therapy needs. You should see your doctor to determine if it is right for you.

EverGo features:

- Weighs less than 10 pounds
- Oxygen capacity 1,050 ml/min
- Non-medical design looks like an ordinary carry-on bag or camcorder case
- Quiet at 43 dBA

Your EverGo system includes:

- **EverGo** portable oxygen concentrator
- Two lithium ion batteries that provide 8 hours of operation*
- AC and DC power cords for charging at home or on the go
- Fold up, lightweight wheeled cart
- Accessories case
- User manual

FAA approved
for air travel



*Based on 20 breaths per minute and a typical device setting of 2



Custom mobile cart
with locking handle



Easy access to batteries
makes them easy to change
and charge – EverGo can
be simultaneously charged
and operated



**Simple-to-use
touch-screen interface**
controls all EverGo functions
and clearly shows battery
life and prescribed settings



Low maintenance
just one easy-to-remove
filter for you to clean

Product specifications

Feature	Specifications
Oxygen concentration	90% +/- 3%
Oxygen delivery pulse dose	Settings 1 - 6 in 0.5 increments
Pulse dose settings	1 - 12 ml, 2 - 23 ml, 3 - 35 ml, 4 - 47 ml, 5 - 59 ml, 6 - 70 ml
Weight	10 lbs (4.5 kg) with 2 batteries, 8.5 lbs (3.9 kg) with 1 battery
Dimensions	Length: 12 in (30.5 cm) Width: 6 in (15.2 cm) Height: 8.5 in (21.6 cm)
Battery duration	2 batteries = 8 hours (setting of 2, 20 bpm) 1 battery = 4 hours (setting of 2, 20 bpm)
Battery recharge time	2 hours – 1 battery 4 hours – 2 batteries (based on fully discharged battery)
Sound level	43 dBA typical (setting of 2, 20 bpm)
Cannula length	30 ft maximum
Operating altitude	Up to 8000 ft
Oxygen purity alarms	Low oxygen at 83%
Alarms and indicators	Low battery at 17 minutes life left; battery status on touch screen and on battery; no breath alarm after 2 minutes; high breath rate alarm
Oxygen capacity	1050 ml per minute (pulse dosed)

Ordering information

Feature	Part number
Concentrator systems	
EverGo system with black carrying case	900-000
EverGo system with navy blue carrying case	900-002
Accessories	
Air inlet filter (6 pack)	900-150
EverGo black carrying case	900-200
EverGo blue carrying case	900-202
Rechargeable Li-ion battery	900-102
Accessory case	900-250
Battery charger/recalibrator	900-103
Mobile cart	900-100
DC power cord	900-101
AC power supply	900-105
AC power cord (American)	900-104
Extended warranty	
EverGo 1-year extended warranty	1038653
EverGo 2-year accident protection	1038570
EverGo 3-year accident protection (includes 1-year extended warranty)	1038648
Maintenance plans	
Deluxe maintenance plan	1038651
Premium maintenance plan	1038652

For more information on EverGo, talk to your homecare provider or visit evergo.respironics.com.



EverGo
portable oxygen
concentrator

**Philips Healthcare is part
of Royal Philips Electronics**

How to reach us

www.philips.com/healthcare

Fax: +31 40 27 64 887

Asia

+852 2821 5888

Europe, Middle East, Africa

+31 40 27 87246

Latin America

+55 11 2125 0764

North America

+1 425 487 7000

800 285 5585

(toll free, US only)

Philips Home

Healthcare Solutions

1010 Murry Ridge Lane

Murrysville, PA 15668

Customer Service

800 345 6443 (toll free)

+1 724 387 4000

Respironics Europe,

Africa, Middle East

+33 1 47 52 30 00

Respironics Asia Pacific

+852 3194 2280

www.philips.com/respironics

Respironics and EverGo are trademarks of Respironics, Inc. and its affiliates.

©2009 Koninklijke Philips Electronics N.V.

All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Geyer MB 03/16/09 MCI 4102100 PN 1059357