



## Designed with the tiniest babies in mind

The right fit as well as comfort

**WeePee** is the only developmentally appropriate diaper for premature infants that comes in three sizes. For babies weighing less than 800 grams to infants weighing 2250 grams, WeePee has you covered. And their soft, easy-to-fasten, non-adhesive tabs ensure a comfortable fit.

**PHILIPS**

Children's  Medical Ventures

## WeePee diapers

WeePee diapers are specifically designed for premature infants to promote increased comfort and developmentally appropriate positioning. Available in three sizes, WeePee's incorporate non-adhesive open/close fasteners and a highly absorbent rectangular pad. The extra small WeePee snugly fits babies weighing less than 800 grams and can absorb up to 25 cubic centimeters of moisture while the refastenable tabs are soft and flexible and won't stick to the baby's skin.

Designed with input from Neonatal Intensive Care Unit (NICU) nurses, and with the tiniest babies in mind, WeePee diapers combine the following features to provide comfort and security.

### Three sizes for a custom, comfortable fit for all preemies

- Extra small (for babies <800 grams)
- Small (for babies 800-1250 grams)
- Medium (for babies 1250-2250 grams)

### Narrow, absorbent pad

- Super absorbent to keep moisture away from baby's skin  
— Amount absorbed: XS: 25cc; S: 50cc; M: 70cc
- Designed to promote developmentally appropriate positioning

### Soft, flexible, refastenable closures

- Easy to attach and reattach diaper tabs for a custom fit
- Will not stick to baby's skin

### Double-stranded, softer elastic leg gathers

- Snugger fit helps minimize leaks
- Softer elastic enhances comfort on delicate skin

### Color-coded packaging to help differentiate sizes at a glance

- Purple for extra small
- Teal for small
- Peach for medium



Extra Small

Small

Medium



Extra small

(for babies <800 grams)

Small

(for babies 800-1250 grams)

Medium

(for babies 1250-2250 grams)

### Ordering information

Item no.	Description	Quantity
1017824	Extra small (XS) <800 g.	360/case*
1017855	Small (S) 800-1250 g.	360/case*
1017856	Medium (M) 1250-2250 g.	360/case*

\*30 diapers per sleeve; 12 sleeves per case

WeePee is a trademark of Respironics, Inc. and its affiliates.



©2009 Koninklijke Philips Electronics N.V. All rights are reserved. Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

GEYER JAH 05.01.09 MCI 4101863 PN 1056110

Respironics Europe,  
Africa, Middle East  
+33 1 47 52 30 00  
Respironics Asia Pacific  
+852 3194 2280

Philips Children's Medical Ventures  
191 Wyngate Drive  
Monroeville, PA 15146  
800 345 6443  
+1 724 387 4000  
[www.philips.com/childmed](http://www.philips.com/childmed)